



FORMULATION

TYGON® F-4040-A

DUROMETER HARDNESS Shore A, 57

Hydrocarbon Resistant TAAT Tubing

- Enables safe handling of most petroleum-based fluids
- Low rate of permeation
- Flexible and easy to install
- Resistant to swelling and embrittlement
- Ozone and UV light resistant

Applications

- Hydrocarbon liquids dispensing
- Sterile filling and dispensing
- Laboratory analytical instruments



2-Stop Autoanalysis Configuration

Saint-Gobain Part Number				I.D.		Wall Thickness		Collar
Tygon S3™ E-3603	Tygon S3™ E-LFL	PharMed® BPT	Tygon® F-4040-A	In.	mm	In.	mm	Color Code
A053401-E	A055501-E	n/a	n/a	0.008	0.19	0.035	0.09	Orange/Red
A053402-E	A055502-E	A053903-E	A053303-E	0.010	0.25	0.035	0.09	Orange/Blue
A053405-E	A055505-E	A053905-E	A053305-E	0.020	0.51	0.035	0.09	Orange/Yellow
A053409-E	A055509-E	A053908-E	A053308-E	0.035	0.89	0.033	0.85	Orange/Orange
A053413-E	A055513-E	A053910-E	A053310-E	0.045	1.14	0.033	0.85	Red/Red
A053416-E	A055516-E	A053912-E	A053312-E	0.056	1.42	0.033	0.85	Yellow/Yellow
A053421-E	A055521-E	A053915-E	A053315-E	0.081	2.06	0.033	0.85	Purple/Purple
A053424-E	A055524-E	A053918-E	A053318-E	0.110	2.79	0.033	0.85	Purple/White

3-Stop Autoanalysis Configuration

Saint-Gobain Part Number				I.D.		Wall Thickness		Collar
Tygon S3™ E-3603	Tygon S3™ E-LFL	PharMed® BPT	Tygon® F-4040-A	In.	mm	In.	mm	Color Code
A053501-E	A055601-E	n/a	n/a	0.008	0.19	0.035	0.09	Orange/Red
A053502-E	A055602-E	A054003-E	A054303-E	0.010	0.25	0.035	0.09	Orange/Blue
A053505-E	A055605-E	A054005-E	A054305-E	0.020	0.51	0.035	0.09	Orange/Yellow
A053509-E	A055609-E	A054008-E	A054308-E	0.035	0.89	0.033	0.85	Orange/Orange
A053513-E	A055613-E	A054010-E	A054310-E	0.045	1.14	0.033	0.85	Red/Red
A053516-E	A055616-E	A054012-E	A054312-E	0.056	1.42	0.033	0.85	Yellow/Yellow
A053521-E	A055621-E	A054015-E	A054315-E	0.081	2.06	0.033	0.85	Purple/Purple
A053524-E	A055624-E	A054018-E	A054318-E	0.110	2.79	0.033	0.85	Purple/White

Nominal Tubing is also available through Saint-Gobain. Contact one of our valued distributors for details.

MANUFACTURING PLANT:

2664 Gilchrist Road, Akron, OH 44305 | Tel: (800) 798-1554 Tel: (330) 798-9240 Fax: (330) 798-6968

Tygon® Autoanalysis Tubing